

## **Certificate in relation to air transportation by passengers with walking disabilities**

We hereby certify that our product

### **e-motion M25 (incl. DuoDrive)**

including the lithium-ion batteries fulfil the requirements of

### **IATA chapter 2.3.2.4**

and are therefore suitable for transportation by passengers according  
The pilot in command must be informed about the location of the mobility aid.

The rated energy of the built-in drive unit battery (which is not intended for removal by  
the user) is **230 Wh**. The rated energy of the DuoDrive control unit is **11 Wh**.

The e-motion M25 can be set into a dedicated "flight mode". In this mode the system  
is completely deactivated. While in this mode the e-motion M25 can only be switched  
on by a coded procedure. Unintended switching on by luggage is not possible.  
Furthermore, the Bluetooth device of the system is completely switched off.  
Since wheels and control unit are portable and contain batteries, we recommend that  
you take them with you to the passenger cabin.

**We strongly recommend contacting the airline sufficiently in advance to  
request approval for transport.**

This certificate does not apply to damaged or defective batteries.

Due to the frequent amendments of the respective regulations, which partly occur on  
an annual basis this certificate is valid until **31<sup>st</sup> December 2024**.

Albstadt, 09<sup>th</sup> of July 2024

Alber GmbH



Joost Beltman  
Managing Director